CELLO™

Balloon Guide Catheter

Intravascular Catheter with Cuff

Instructions for Use

Balloon Guide Catheter	Instructions for Use	EN	4
Cathéter-guide à ballonnet	Mode d'emploi	FR	6
Ballon-Führungskatheter	Gebrauchsanweisung	DE	8
Catetere guida a palloncino	Istruzioni per l'uso	IT	10
Catéter guía con globo	Instrucciones de uso	ES	12
Ballonguidekateter	Bruksanvisning	SV	14
Ballon Gids Katheter	Gebruiksaanwijzing	NL	16
Cateter-guia com Balão	Instruções de utilização	PT	18
Ballonki-katetri	Käyttöohjeet	FI	20
Ballonguidekateter	Brugsanvisning	DA	22
Οδηγός Καθετήρας με Μπαλόνι	Οδηγίες χρήσης	EL	24
Balónkový vodící katétr	Návod k použití	CS	26
Ballonos vezető katéter	Használati előírások	HU	28
Баллонный направляющий катетер	Инструкция по применению	RU	30
Balonowy cewnik prowadzący	Instrukcja użytkowania	PL	32
Balon Kılavuz Kateter	Kullanım Talimatı	TR	34
Ballongledekater	Bruksanvisning	NO	36
Balónkový vodiaci katéter	Návod na použitie	SK	38
Cateter cu ghidaj cu balon	Instrucțiuni de utilizare	RO	40
Водещ балонен катетър	Инструкции за употреба	BG	42
Ballooniga juhtekateeter	Kasutusjuhised	ET	44
Balonski uvodni kateter	Upute za upotrebu	HR	46
Balioninis kreipiamasis kateteris	Naudojimo instrukcija	LT	48
Vadīšanai paredzēts balonkatetrs	Lietošanas instrukcija	LV	50
Balonski vodilni kateter	Navodila za uporabo	SL	52
Cateter-guia com Balão	Instruções de Utilização	PT- BR	54
Ống thông can thiệp đầu gắn bóng	Hướng dẫn sử dụng	VI	56
Балонски уводни катетер	Упутство за употребу	SR	58
Шарлы бағыттағыш катетер	Пайдалану жөніндегі нұсқаулық	KK	60
Балон катетер за жица водич	Упатство за употреба	MK	62
Kateter Pemandu Balon	Petunjuk Penggunaan	ID	64

Device Description

The Cello™ Balloon Guide Catheter is a coaxial-lumen, braid-reinforced, variable stiffness catheter with two radiopaque markers on both the distal and proximal ends of the balloon and a bifurcated luer hub on the proximal end. A compliant silicone balloon is mounted on the distal end. Balloon Guide Catheter dimensions and recommended balloon inflation volumes are indicated on the product label. If indicated on product label, a dilator is provided.

Indications for Use

The CELLO Balloon Guide Catheter is indicated for use in facilitating the insertion and guidance of intravascular catheters into a selected blood vessel in the peripheral and neuro vasculature systems. The balloon provides temporary vascular occlusion during these and other angiographic procedures.

Contraindications

Use of the Cello™ Balloon Guide Catheter is contraindicated under these circumstances

- Patients having acute cardiac infarction, severe irregular heartbeat or abnormal serum electrolyte.
- Patients having severe disorder of liver or kidney function.
- Patients having a blood coagulation disorder or severe changes in coagulation for any reason.
- Patients who cannot lie supine on the angiography bed because of congestive heart failure or respiratory problems.
- Patients having irritation against contrast media.
- Pregnant women.
- Patients who are uncooperative during examination.

Complications

Possible complications include, but are not limited to, the following: infection, hematoma, distal embolization, vessel thrombosis, dissection, false aneurysm formation, acute occlusion, clot formation, hemorrhage at the puncture site, intracranial hemorrhage, arterial rupture, arterial injury, arterial perforation, Arteriovenous malformation, infection and pain, nausea and vomiting, blood pressure drop, arrhythmia, late bleeding, vasospasm, obstruction of a blood vessel, nervous system injury, stroke and death.

Compatibility

Introducer sheath French size must be greater than or equal to balloon guide catheter French size. Maximum guidewire diameter is indicated on the product label.

Warnings

- The CELLO Balloon Guide Catheter should only be used by physicians who have received appropriate training in interventional techniques.
- Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning.
- Never advance or torque catheter against resistance without careful assessment of cause of resistance using fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in damage to vessel or catheter.
- To reduce risk of complications due to slow balloon deflation, adhere to the following recommendations:
 - Wet the distal shaft with saline before passing into the introducer sheath.
 - Minimize pushing forces on the shaft during advancement. These forces can cause wrinkles in the shaft that can slow balloon deflation.
 - Do not use the device if the shaft is damaged during use.
 - Prepare balloon according to the recommended procedure.
- To reduce the risk of complications due to air emboli, remove air from the balloon according to the recommended procedure.
- Withdrawing the balloon through the introducer sheath may damage the balloon. Do not use the catheter again after withdrawing the balloon through the introducer sheath.
- To avoid balloon leakage, do not allow the balloon to contact calcified or stented arteries and do not allow the balloon to move during inflation.
- Do not use a device that has been damaged. Use of damaged devices may result in complications.
- Do not exceed maximum recommended balloon inflation volume as indicated on the label. Excess inflation volume may rupture the balloon.
- For through-lumen, do not exceed 2.07 MPa (300 psi) maximum recommended infusion pressure. Excess pressure may result in catheter rupture or tip detachment.
- If flow through the catheter becomes restricted do not attempt to clear catheter lumen by infusion. Doing so may cause the catheter to rupture, resulting in vessel trauma. Remove and replace the catheter.
- Do not steam shape the guide catheter.
- Do not pull the catheter at a force exceeding 3.5 N (0.35 kgf).
- Balloon should be inflated by one of the followings;
 1. Inflating the balloon with the dilator inserted (recommended)
 2. Inflating the balloon 3 seconds or more with 1 mL syringe.

Precautions

- Single use only. (If the product is to be re-used, it may cause infection to patients or damage to the product.)
- Store in a cool, dry, dark place.
- Do not use open or damaged packages.
- Use by “Use By” date.
- Exposure to temperatures above 54°C (130°F) may damage the device and accessories. Do not autoclave.
- Upon removal from package, inspect the device to ensure it is not damaged.
- Do not expose the device to solvents.
- Do not use oily contrast media, consisting of ethyl ester of iodinated poppy-seed oil fatty acids, such as Lipiodol.
- The guidewire used with this system should incorporate a hydrophilic coating.
- Use the device in conjunction with fluoroscopic visualization and proper anticoagulation agents.
- Do not conduct MRI inspection while the product is in place.
- Torquing the guide catheter while kinked may cause damage that could result in separation of the catheter shaft.
- If a device becomes lodged in the guide catheter, or if guide catheter becomes severely kinked, withdraw the entire system (guide catheter, guidewire and catheter sheath introducer).
- To prevent thrombus formation and contrast media crystal formation maintain a constant infusion of appropriate flush solution through the guide catheter lumen.
- After use, dispose of this product and packaging based on the hospital, administrative and /or local government policy.
- If user and/or patient notice any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority.

Recommended Procedure

1. Gently remove the catheter from the pouch by grasping the hub and slowly pulling the catheter straight out of the protective tubing without bending the shaft. Inspect the catheter before use to verify that its size, shape and condition are suitable for the specific procedure.
2. Dilator Instructions (if applicable):
 - Remove dilator from protective tube.
 - Gently insert dilator into guide catheter lumen.
 - Flush dilator with heparinized saline.
3. Prepare balloon inflation media by mixing contrast with saline (50% by volume or equivalent to 150 mg/mL iodine).
4. Fill a 20 mL syringe with 5 mL of balloon inflation media. Attach a 3-way stopcock to the balloon hub. Attach a 20 mL syringe to the stopcock.
5. Attach a 1 mL syringe to the balloon’s 3-way stopcock. Turn the stopcock towards the 1 mL syringe.
6. With the 20 mL syringe pointing downward:
 - Pull back on the syringe plunger to aspirate the balloon lumen. Maintain negative pressure until the air bubbles stop forming in the syringe.
 - Release the syringe plunger to allow media to be drawn into the balloon lumen. Do not infuse media.
 - Again pull back on the syringe plunger to aspirate the balloon lumen. Maintain negative pressure until air bubbles stop forming in the syringe.
 - Release the syringe plunger to allow media to be drawn into the balloon lumen. Do not infuse media.
7. Turn stopcock towards the balloon hub. Transfer maximum recommended balloon inflation volume from the 20 mL syringe to the 1 mL syringe.
8. Inflate the balloon up to the recommended inflation volume. Refer to the product label. Turn the stopcock off towards the balloon hub.
9. Inspect the balloon for leakage. Keep the balloon inflated until the air bubbles diffuse from the balloon. If air remains in the balloon, it should dissipate if left in a sterile field for several minutes.
10. Deflate the balloon by turning stopcock towards the 1 mL syringe and aspirating with the 20 mL syringe.
11. After ensuring that the balloon is fully deflated, wet the distal shaft with saline and insert the balloon portion of the catheter into the inserter.
12. Gently insert the tip of the guidewire and the inserter/guide catheter assembly through the proximal valve of the introducer sheath.
13. During insertion of the catheter, when the balloon passes the sheath inducer, pull back the inserter only and peel it away from the catheter.
14. Place the guide catheter in the selected vessel using fluoroscopy.
15. Remove dilator (if applicable) and guidewire.
16. Attach a rotating hemostatic valve (RHV) to the guide catheter through-lumen hub. Flush the through-lumen with heparinized saline. Attach a 3-way stopcock to the RHV and to appropriate flush solution.
17. To inflate the balloon, transfer up to the maximum recommended balloon inflation volume from the 20 mL syringe to the 1 mL syringe and gently infuse the balloon inflation media with the 1 mL syringe until desired balloon diameter is attained.
18. Ensure the balloon is completely deflated before withdrawing the guide catheter.

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MicroTherapeutics, Inc.
Fuji Systems Corporation Shirakawa Plant
MedDevConsult GmbH