ENDOTRACHEAL TUBE UNIVENT™

3.5mm, without cuff

4.5mm, with cuff

Instructions for Use

ENDOTRACHEAL TUBE		
UNIVENT™, 3.5mm, without cuff	EN	4
UNIVENT™, 4.5mm, with cuff		
Instructions for use		
SONDE ENDOTRACHÉALE		
UNIVENT™, 3,5 mm, sans ballonnet	FR	7
UNIVENT™, 4,5 mm, avec ballonnet		
Mode d'emploi		
ENDOTRACHEALTUBUS		
UNIVENT™, 3,5 mm, ohne Manschette	DE	10
UNIVENT™, 4,5 mm, mit Manschette		
Gebrauchsanweisung		
TUBO ENDOTRACHEALE		
UNIVENT™, 3,5 mm, senza cuffia	IT	13
UNIVENT™, 4,5 mm, con cuffia		
Istruzioni per l'uso		
TUBO ENDOTRAQUEAL		
UNIVENT™, 3,5 mm, sin manguito	ES	16
UNIVENT™, 4,5 mm, con manguito		
Instrucciones de uso		
ENDOTRAKEAL TUB		
UNIVENT™, 3,5 mm, utan kuff	SV	19
UNIVENT™, 4,5 mm, med kuff		
Bruksanvisningar		
ENDOTRACHEAL TUBE		
UNIVENT™, 3,5 mm, zonder cuff	NL	22
UNIVENT™, 4,5 mm, met cuff		
Gebruiksaanwijzing		
TUBO ENDOTRAQUEAL		
UNIVENT™, 3,5 mm, sem cuff	PT	25
UNIVENT™, 4,5 mm, com cuff		
Instruções de utilização		
ENDOTRAKEALT RØR		
UNIVENT™, 3,5 mm, uden manchet	DA	28
UNIVENT™, 4,5 mm, med manchet		
Brugsanvisning		
ΕΝΔΟΤΡΑΧΕΙΑΚΟΣ ΣΩΛΗΝΑΣ		
UNIVENT™, 3,5mm, χωρίς cuff	EL	31
UNIVENT™, 4,5mm, με cuff		
Οδηγίες χρήσης		
ENDOTRACHEÁLNÍ TRUBICE		
UNIVENT™, 3,5 mm, bez manžety	CS	34
UNIVENT™, 4,5 mm, s manžetou		
Návod k použití		
RURKA DOTCHAWICZA		
UNIVENT™, 3,5mm, bez mankietu	PL	37
UNIVENT™, 4,5mm, z mankietem		
Instrukcja użycia		
ENDOTRAKEAL TÜP		
UNIVENT™, 3,5 mm, kafsız	TR	40
UNIVENT™, 4,5 mm, kaflı		
Kullanım Talimatı		
ENDOTRAKEAL SLANGE		
UNIVENT™, 3,5 mm, uten mansjett	NO	43
UNIVENT™, 4,5 mm, med mansjett		
Bruksinstruksjoner		
ENDOTRAHEALNA CIJEV		
UNIVENT™, 3,5 mm, bez balončica	HR	46
UNIVENT™, 4,5 mm, s balončićem		
Upute za uporabu		

Indications for Use :

UNIVENT™ (TCB Type) is indicated for use in airway management of surgical patients to perform one-lung ventilation.

Warnings / Precautions :

[General]

- Single use only (If the product is to be re-used, it may cause infection to patients or damage to the product).
- This product contains phthalic acid in the Blocker pilot balloon, but it does not contact with the body.
- It should not be replaced with a new piece of the same type after 30 days of use.
- Do not use if the sterile package is open, wet or damaged.
- The tube is for oral intubation only. Do not use for nasal intubation.
- Prior to use, take off and discard the aeration plate mounted at the one way valve on the blocker tube. This aeration plate is used for ventilation to avoid the risk of inflation of the blocker cuff by negative pressure during sterilization. If the aeration plate is detached from one way valve before opening the package, check if the blocker cuff is functioning.
- When pushing blocker cuff out from main channel or drawing back blocker shaft into the main channel, it is suggested to manipulate blocker shaft close to the blocker mantle tube, because 1) blocker cuff tends to be caught on the inner wall at the tip part of main tube and 2) blocker shaft may kink due to lack of stiffness.
- Avoid contact with laser beam or an electrosurgical knife in the immediate area of UNIVENT™. Such contact can result a sudden ignition of UNIVENT™ in the presence of mixtures of nitrous oxide and oxygen or pure oxygen.
- Care must be taken to avoid occlusion at the tip part of main tube and blocker tube when applying lubricant.
- Respiratory circuits of the respirator must be fit with standard 15mm connector. The 15mm connector of this product must be fit firmly with the respiratory circuits.
- Do not cut the tube to adjust its length.
- Unlike the ordinary UNIVENT™, blocker tip of 3.5mm & 4.5mm UNIVENT™ are closed end. Lung collapse and suction cannot be done through the blocker lumen. Collapse the desired lung naturally. Do not disconnect the tube from ventilator during the procedure.
- After use, dispose of this product and packaging based on the hospital, administrative and /or local government policy.
- If user and/or patient notice any serious incident that has occurred in relation to the device, please report to the manufacturer and the competent authority.

[Cuff Related]

- Test the cuff, and valve assembly prior to use. If dysfunction is detected in any part of the inflation system, the tube should not be used.
- After checking blocker cuff, deflate blocker cuff and pull blocker shaft back to the setting end mark (a double line marked on blocker shaft) and have blocker cuff placed at the tip part of the main tube. This is to avoid the bent tip of blocker cuff to narrow the inner lumen of main tube and it also prevents irritation to the lateral wall of trachea and bronchus during the insertion.
- Do not overinflate the blocker cuff. Overinflation can result in bronchus damage. Inflate the blocker cuff by using Luer-Slip type syringe.
- The cuff wall should not be damaged during the insertion of the tube. If any damage is detected, the tube should not be used.
- Apply lubricant on blocker cuff in order to reduce the risk of damage of blocker cuff caused by the friction with the pocket part of the main tube.
- Iodic disinfectant should not be used. It may deteriorate the material of the cuff.
- Inflate the cuff until suitable seal can be reached.
 - Maximum injected air volume for main cuff:
 - I.D. size 4.5mm (with cuff) max. 5mL
 - Note: I.D. size 3.5mm (without cuff) does not have main cuff.
 - Maximum Injected air volume for blocker cuff per size:
 - max. 2mL
- Cuff inflation and pressure should be monitored at all time. Permeation of mixtures of nitrous oxide (laughing gas) and oxygen can cause the change of cuff volume and increase or decrease pressure of silicone cuff.
- Air in the cuff must be removed completely before adjusting the position of main tube cuff or blocker cuff. It may damage the trachea, bronchus, or cuff itself while manipulating without completely deflating of the cuff.
- Before retracting the blocker cuff into the main channel, make sure that the air has been removed from the blocker cuff completely. Then gently pull back the blocker shaft until the setting end mark (a double-line mark) on blocker appears at the end of blocker mantle tube. Do not pull the blocker tube with excessive force. It may damage the blocker cuff.

[Clinical Related]

1. Always use a flexible bronchofiberscope to position and set the tube blocker.
2. The bronchus should be blocked only after the patient has been positioned on the side, the thoracotomy started and the pleura opened so as to reduce the amount of pulmonary blood flow to the upper lung.
3. The blocker cuff pressure should be checked with cuff pressure gauge to prevent the overinflation of the cuff.
4. When the bronchus is blocked the inspired oxygen concentration should be increased to 50% or more, and the patient mechanically ventilated.
5. PaO₂ should be measured when the blocked lung is either visually confirmed to be collapsed, or twenty minutes have elapsed since the blocking of the bronchus. The value of PaO₂ should not decrease thereafter owing to "Hypoxic Pulmonary Vasoconstriction". However it is desirable to check the PaO₂ again to make sure that one lung anesthesia is correctly taking place. In addition, SaO₂ should be monitored with a pulse oximeter.

Recommended Procedure :

1. Remove the sterile UNIVENT™ from the package with care.
2. Prior to use, take off and discard the aeration plate mounted at one way valve on blocker tube.
3. Test inflation of the main cuff and blocker cuff prior to use.
DO NOT USE Luer-Lock syringe; must use Luer-Slip syringe to inflate the cuff.

Maximum injected air volume for main cuff:

I.D. size 4.5mm (with cuff) max. 5mL

Note: I.D. size 3.5mm (without cuff) does not have main cuff.

Maximum Injected air volume for blocker cuff per size:
max. 2mL

4. Deflate air in main cuff and blocker cuff completely after the inflation test.
5. Apply lubricant on both cuffs.
6. Pull the blocker shaft back slowly to the setting end mark (a double-line marked on the blocker shaft) and have blocker cuff placed inside the main channel.
7. Move the cap stopper while give it a twist by finger and put it into the end of the blocker mantle tube.
8. Fix the blocker with the band stopper on main tube.
9. When connecting to respirator, make sure the 15mm connector is firmly attached to the respiratory circuit.
10. Intubate UNIVENT™ by referring the following methods:

[Main Tube Rotation Method]

- A. UNIVENT™ is inserted into the trachea.
- B. Turn the tube 90 degrees towards the operative side, this may be facilitated by turning the patient's head; to check the degree to tube rotation make sure the tube exits the patient's mouth on the side for thoracotomy. Rotate the tube with the tracheal cuff deflated. After rotation, inflate trachea cuff and fix the tube firmly at patient mouth with cloth tape.
- C. Push the blocker shaft; the blocker will follow the lateral wall of the trachea and enter the target main stem bronchus. Use a flexible bronchofiberscope to position blocker as deeply as possible to prevent cuff herniation into the carina and to prevent blocker dislodgment when the patient is turned on their side. Following this, fix the blocker to the endotracheal tube with the cap stopper and band stopper.
- D. After the patient has been turned on the side and thoracotomy has taken place, place 5mL of air into the blocker cuff, checking cuff pressure. Pitting of the lung surface which appears and remains following finger compression is evidence that the bronchus is completely blocked.

[Blocker Rotation Method]

- A. Insert UNIVENT™ and have the tube fixed.
 - B. Insert bronchofiberscope into the tube and then push while twisting blocker tube to the target main stem bronchus under direct vision. The rest of the manipulations are same as Main Tube Rotation Method.
11. When extubating UNIVENT™, make sure that the air has been removed from the blocker cuff, then gently pull back the blocker until the setting end mark (a double-line mark) on blocker appears at the end of blocker mantle tube and the blocker cuff is placed in the pocket at the tip part of the main tube. Then deflate the main cuff completely. Extubation should be performed currently accepted medical techniques.

Adverse Reactions :

The following adverse reactions have been reported to be associated with the use of UNIVENT™ during the intubation procedure, the intubation period, or in extubation procedure. The order of listing is alphabetical and does not indicate frequency or severity: cartilage necrosis; consequences of failure to ventilate including death; damage to the perichondrium; emphysema; excoriated membranes of pharynx; glottic edema; infections; laryngeal obstruction; laryngeal stenosis; bronchitis; submucosal hemorrhage; tracheorrhagia; tracheal stenosis; traumas (lips, pharynx, trachea, glottis, and etc.).

Storage Conditions :

Keep the product dry and store in clean conditions, avoiding high temperature, humidity and direct sunlight.

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MedDevConsult GmbH
HSC-Medical GmbH
Fuji Systems Corporation Shirakawa Plant